


Ponto™ – Bone Anchored Hearing System

Ponto 5 Mini More sound to the brain



oticon
MEDICAL

A woman with curly hair is hugging a young child from behind on a city street. The woman has her eyes closed and a peaceful expression. The child is wearing a blue shirt and is also hugging the woman. The background shows a modern city street with buildings and a clear sky. There are some faint, glowing lines and circles overlaid on the image, suggesting a digital or technological theme.

**Enjoy the world's
smallest sound processor
– it's still a Ponto**

Whether you are new to Ponto or want to upgrade your existing sound processor, Ponto 5 Mini gives you access to more sound and proven better speech understanding with the OpenSound Navigator.^{1,2}

Packed with technology

- Small and discreet
- Enjoy the full soundscape
Hear sounds all around you with OpenSound Navigator™ technology
- Premium sound quality
With no audible feedback due to the patented OpenSound Optimizer™*1
- Listening comfort
Tailored programmes for different situations
- Part of your digital life
Simple to connect to your favourite devices
- Remote appointments
Get your device adjusted online



*When fitted according to recommendation

Stay connected with your sound environment

Conventional sound processors aim to improve speech understanding by removing sound, for example by using directionality to just focus on the speakers in front of you. While this may improve speech understanding, it can also leave you feeling cut off.

Ponto 5 Mini takes a different approach. Its pioneering OpenSound Navigator technology lets you hear sounds from all around you.



This makes it easy to follow what's going on. And the results of this technology speak for themselves: speech understanding is improved by more than 20%.²

In addition to creating an open sound environment, the pioneering Ponto 5 Mini also gives you more stable sound and prevents audible feedback.¹

If you're experiencing feedback with your current bone anchored system, ask your hearing care professional to try a Ponto 5 Mini.






Try the open sound experience for yourself.
Ask for a trial of Ponto 5 Mini.

Better speech understanding with less listening effort


Ponto 5 Mini uses OpenSound Navigator technology to open up the full soundscape so you can hear sounds from all directions.

New research confirms that this open soundscape gives the brain what it needs to create a better listening experience. It has been proven to give better speech understanding with less listening effort.²

Only Ponto bone anchored sound processors with BrainHearing™ technology offering open sound experience give you access to the full soundscape.



Connect, stream and relax



Whether you want to have a video call, watch TV or talk on the phone, your Ponto 5 Mini can wirelessly connect to a wide range of digital devices used in daily life.

- **2.4 GHz Bluetooth Low Energy**
connect to your smartphone
- **Oticon ON App**
adjust the volume, change programmes and much more using your phone
- **EduMic**
stream the teacher's voice or connect to FM or telecoil



- **ConnectClip**
stream to and from your smartphone and use as a remote microphone
- **TV Adapter 3.0**
stream sound from your TV
- **Oticon RemoteCare App**
attend appointments with your audiologist online
- **IFTTT – the Internet of Things**
connect to internet-enabled devices and services
- **Ponto Care App**
support to get the best experience with Ponto 5 Mini



The ultimate in discretion and reliability

Ponto 5 Mini is small, discreet and available in different colours to blend in with your hair. Or you can add your personal style using the bright and colourful skins or stickers.

Performance you can rely on

All aspects of Ponto 5 Mini have undergone extensive testing, including drop tests, humidity and sweat tests, extreme temperature changes and full immersion in water.* That's how you can be sure that the Ponto 5 Mini won't let you down.

LED indicator

The LED indicator light can be activated to help caregivers monitor the status of the processor.



I was incredibly surprised two weeks later when I had the final fittings. That I was able to hear so much and from all around me. I was amazed.”

Fatima, Ponto user, France

*IP 57 classification

References:

¹Data on file at Oticon Medical

²Manuscript in preparation, Data on File - Clinical study BC102

Because sound matters

Oticon Medical is a global company in implantable hearing solutions, dedicated to bringing the power of sound to people at every stage of life. As part of the Demant group, a global leader in hearing healthcare with more than 16,500 people in over 30 countries and users benefitting from our products and solutions in more than 130 countries, we have access to one of the world's strongest research and development teams, the latest technological advances and insights into hearing care.

Our competencies span more than a century of innovations in sound processing and decades of pioneering experience in hearing implant technology. We work collaboratively with patients, physicians and hearing care professionals to ensure that every solution we create is designed with users' needs in mind. We have a strong passion to provide innovative solutions and support that enhance quality of life and help people live full lives – now and in the future.

Because we know how much sound matters.



 **Oticon Medical AB**
Datavägen 37B
SE-436 32 Askim
Sweden
Tel: +46 31 748 61 00

231708UK / 2021.06