

Oticon Medical User Order Form

CONTACT INFORMATION	Today's Date _____
Contact Name _____	
Contact Email _____	
Contact Phone (___) _____ - _____	Send Email Confirmation Yes No





SHIP TO	<i>(Please complete all information including name and phone number)</i>
All orders ship via \$10 flat rate UPS Next Day.	
Name _____	
Address (line 1) _____	
Address (line 2) _____	
City _____	State _____ Zip _____
Phone (___) _____ - _____	

PATIENT INFORMATION	
First Name _____	<i>Middle Initial</i> _____
Last Name _____	
Date of Birth ____ / ____ / ____ <small>MM DD YYYY</small>	





Products

Wireless Accessories









The below accessories are compatible with the Ponto 5 Mini, Ponto 5 SuperPower, Sentio 1 Mini, and Ponto 4 sound processors.

Product Image	Item Number	Description	Price	Quantity
	178509 ConnectClip	ConnectClip combines a Bluetooth® headset and a microphone in one. ConnectClip streams sound directly from a wireless device, e.g. a smartphone to the sound processor in high quality while also functioning as a wireless microphone and remote control.	\$220	_____
	169229 Remote Control 3.0	Remote Control 3.0 offers discreet and easy control over your sound processor. Use it to adjust volume, change programs or mute your devices.	\$160	_____
	168628 TV Adapter 3.0	TV Adapter 3.0 connects your sound processors with the TV for wireless streaming of the TV sound directly to the sound processor.	\$170	_____
	191009 EduMic	EduMic is an educational solution which allows teachers to stream voice directly to several hearing device users in a classroom. It also can connect to most devices that feature an AC jack and provides access to telecoil mode.	\$725	_____










Care & Maintenance

Product Image	Item Number	Description	Price	Quantity
	173365 Safety line	The safety line may be attached to the your Ponto or Sentio and the patient's clothing to reduce risk of loss or damage to the sound processor should it come loose.	\$5.40	_____
	890-60-232-00 MultiTool	The back-end of the MultiTool is magnetic and can be used for removing and changing batteries.	\$1.09	_____
	19900-4020 Mini Super Dri-Aid Kit	Hearing aid dehumidifier jar for removing moisture from sound processors.	\$12	_____
	19900-4080 Zephyr Dry and Store Kit	Electronic dehumidifier for your sound processors. Uses heat, moving air, and a desiccant to extract moisture from your device.	\$89	_____






Sentio Specific Accessories

Product Image	Item Number	Description	Price	Quantity
	265231 Sentio Side Indicator	The Sentio side indicator is a plug kit containing plugs for side indication in the Sentio 1 Mini sound processor. Each kit contains two blue and two red plugs.	\$4	_____
	230821 TAR Tool	The TAR tool is used for opening and closing the tamper-resistant lock on the sound processor and to lock or unlock the magnet on Sentio 1 Mini sound processors.	\$0	_____
	265232 TAR Battery lid set, Chroma Beige	Each set includes three lids with TAR (tamper-resistant) hole for locking.	\$51.32	_____
	265232 TAR Battery lid set, Terracotta		\$51.32	_____
	265234 TAR Battery lid set, Silver		\$51.32	_____
	265235 TAR Battery lid set, Steel Grey		\$51.32	_____
	265236 TAR Battery lid set, Chestnut Brown		\$51.32	_____
	265237 TAR Battery lid set, Diamond Black		\$51.32	_____






Ponto Specific Accessories

Product Image	Item Number	Description	Price	Quantity
	M50294 Abutment cover, black	The abutment cover can be used to cover the abutment when the sound processor is not connected.	\$2.50	_____
	M50290 Abutment cover, beige		\$2.50	_____
	M51177 Ponto Care kit	The Ponto Care kit can be used to facilitate cleaning of the Ponto System and to establish a good daily aftercare routine. This kit includes wet wipes, cleaning brush, Q-tips, abutment cover, mirror, a mini super dry aid kit, and an aftercare guide.	\$0	_____
	19900-4000 Pediatric After Care Kit	The Pediatric Ponto Care Kit can be used to facilitate cleaning of the Ponto System and to establish a good daily aftercare routine. This kit includes a care kit with wet wipes, cleaning brush, Q-tips, abutment cover, and mirror, a mini super dry aid kit, an aftercare guide, a plush toy, and a Ponto canvas carrier bag.	\$0	_____
	M50030 Cleaning brush	The soft cleaning brush is used to clean the outside and inside of the abutment.	\$1	_____
	204265 Plug kit, Ponto 5 Mini & Ponto 4	The plug kit contains plugs for covering the programming socket on the sound processor and for side indication. Each kit contains three blue, three red and three black plugs.	\$11	_____
	226490 Plug kit, Ponto 5 SuperPower		\$11	_____
	204269 Battery drawer kit, Nail grip for Ponto 5 Mini & Ponto 4	The kits for Ponto 4 and Ponto 5 Mini come in either a nailgrip or in a tamper-resistant version. The battery drawers for Ponto 5 SuperPower only come in one version with a built-in tamper-resistant mechanism. Each kit contains three battery drawers and a tool for removing the drawer.	\$32	_____
	204268 Battery drawer kit, Tamper-Resistant, Ponto 5 Mini & Ponto 4		\$32	_____
	226557 Battery drawer kit, Tamper-Resistant, Ponto 5 Super Power		\$32	_____
	165593 TAR Tool	A tool for opening the tamper-resistant battery drawer on Ponto 4 and Ponto 5 Mini.	\$0	_____
	M50302 Test Rod	The test rod is used for listening, testing and demonstrating the sound processor without connecting it to an abutment. Simply connect the sound processor to the test rod and press it against the head.	\$0	_____






Stickers & Skins for the Ponto 4 & 5 Mini

Product Image	Item Number	Description	Price	Quantity
	197957 Stickers	The adhesive stickers can easily be applied to the sound processor to customize its appearance.	\$6	_____
	234062 Skins, Turquoise/Pink	The skin can be easily placed on the sound processor to change the appearance. Each set consists of six skins in two colors, three of each color. The skins for Ponto Mini fit both Ponto 4 and Ponto 5 Mini.	\$29	_____
	206535 Skins, Red/Yellow		\$29	_____
	234063 Skins, Light Green/Navy		\$29	_____
	234064 Skins, Light Blue/Dark Purple		\$29	_____

Stickers & Skins for the 5 SuperPower

Product Image	Item Number	Description	Price	Quantity
	197957 Stickers	The adhesive stickers can easily be applied to the sound processor to customize its appearance.	\$6	_____
	225096 Turquoise/Pink	The skin can easily be placed on the sound processor to change the appearance. Each set consists of six skins in two colors, three of each color.	\$29	_____
	225098 Skins, Light Green/Navy		\$29	_____
	225099 Skins, Light Blue/Dark Purple		\$29	_____
	225097 Skins, Red/Yellow		\$29	_____

Stickers & Skins for the Ponto 3 Family

Product Image	Item Number	Description	Price	Quantity
	M52857	The adhesive stickers can easily be applied to the sound processor to customize its appearance.	\$5	_____
	Stickers, Ponto 3, Right			
	M52858		\$5	_____
	Stickers for Ponto 3, Left			
	M52859		\$5	_____
	Stickers for Ponto 3 Power & SuperPower, Right			
	M52860		\$5	_____
	Sticker for Ponto 3 Power & SuperPower, Left			
	M52561	The skin can be easily placed on the sound processor to change the appearance. Each set consists of five skins in different colors: cerise, turquoise, orange, navy blue and yellow.	\$29	_____
	Skins, Ponto 3, Right			
	M52562		\$29	_____
	Skins, Ponto 3, Left			
	M51844	The skin can be easily placed on the sound processor to change the appearance. Each set consists of five skins in different colors: pink, green, red, light blue and purple.	\$29	_____
	Skins, Ponto 3, Right			
	M51845		\$29	_____
	Skins, Ponto 3, Left			
	M52563	The skin can be easily placed on the sound processor to change the appearance. Each set consists of five skins in different colors: cerise, turquoise, orange, navy blue and yellow.	\$29	_____
	Skins, 3 Power & SuperPower, Right			
	M52564		\$29	_____
	Skins, 3 Power & SuperPower, Left			
	M51846	The skin can be easily placed on the sound processor to change the appearance. Each set consists of five skins in different colors: pink, green, red, light blue and purple.	\$29	_____
	Skins, 3 Power & SuperPower, Right			
	M51847		\$29	_____
	Skins, 3 Power & SuperPower, Left			

Softbands & Headbands

The softband and headband devices are intended for users who have not had implant surgery but can benefit from a bone conduction solution. The softband consists of an adjustable elastic band with one or two Connectors for the sound processor(s). It has a safety release feature, designed to open if the band gets caught. The softbands are available for unilateral and bilateral use.








Softband 5 and Connector Pad

The Softband 5 comes in two sizes – XS/S and M/L. Each softband pack contains Connector Pads – three Connector Pads for the unilateral band and six Connector Pads for the bilateral bands.

All unilateral softbands are \$65 and all bilateral softbands are \$86.









Softband 5, unilateral XS/S for head
circumference: 380-500mm



		Qty.
	227690 Softband 5, unilateral XS/S, black	_____
	227691 Softband 5, unilateral XS/S, dark brown	_____
	227692 Softband 5, unilateral XS/S, light brown	_____
	227693 Softband 5, unilateral XS/S, beige	_____
	227694 Softband 5, unilateral XS/S, light blue	_____
	227695 Softband 5, unilateral XS/S, pink	_____
	227696 Softband 5, unilateral XS/S, navy blue	_____
	227697 Softband 5, unilateral XS/S, red	_____

Softband 5, unilateral M/L for head
circumference: 500-600mm



		Qty.
	227686 Softband 5, unilateral M/L, black	_____
	227687 Softband 5, unilateral M/L, dark brown	_____
	227688 Softband 5, unilateral M/L, light brown	_____
	227689 Softband 5, unilateral M/L, beige	_____
	267634 Softband 5, unilateral M/L, light blue	_____
	267635 Softband 5, unilateral M/L, pink	_____
	267636 Softband 5, unilateral M/L, navy blue	_____
	267637 Softband 5, unilateral M/L, red	_____

Softband 5, bilateral XS/S for head
circumference: 380-500mm

		Qty.
	227702 Softband 5, bilateral XS/S, black	_____
	227703 Softband 5, bilateral XS/S, dark brown	_____
	227704 Softband 5, bilateral XS/S, light brown	_____
	227705 Softband 5, bilateral XS/S, beige	_____
	227706 Softband 5, bilateral XS/S, light blue	_____
	227707 Softband 5, bilateral XS/S, pink	_____
	227708 Softband 5, bilateral XS/S, navy blue	_____
	227709 Softband 5, bilateral XS/S, red	_____

Softband 5, bilateral M/L for head
circumference: 500-600mm

		Qty.
	227698 Softband 5, bilateral M/L, black	_____
	227699 Softband 5, bilateral M/L, dark brown	_____
	227700 Softband 5, bilateral M/L, light brown	_____
	227701 Softband 5, bilateral M/L, beige	_____
	267638 Softband 5, bilateral M/L, light blue	_____
	267639 Softband 5, bilateral M/L, pink	_____
	267640 Softband 5, bilateral M/L, navy blue	_____
	267641 Softband 5, bilateral M/L, red	_____



226779
Connector Pad set
(pack of 20 pcs)

Connector Pads are used with Softband 5 for wearing comfort and to optimize fit and transmission. The Connector Pads may be replaced as needed.

\$32

All unilateral softbands are \$65 and all bilateral softbands are \$86.

Ponto Softband and Ponto SoundConnector™















The Ponto Softband has a narrower width (15mm) compared to Softband 5. It is available in 14 colors and for unilateral or bilateral use.

Ponto Softband, unilateral



Ponto Softband, bilateral



	Qty:
 173366 Ponto Softband, beige	_____
 173367 Ponto Softband, yellow	_____
 173368 Ponto Softband, red	_____
 173369 Ponto Softband, pink	_____
 173379 Ponto Softband, cerise	_____
 173380 Ponto Softband, purple	_____
 173381 Ponto Softband, navy	_____
 173382 Ponto Softband, indigo	_____
 173383 Ponto Softband, blue	_____
 173384 Ponto Softband, turquoise	_____
 173385 Ponto Softband, green	_____
 173386 Ponto Softband, grey	_____
 173387 Ponto Softband, brown	_____
 173388 Ponto Softband, black	_____

	Qty:
 185311 Ponto Softband, bilateral, beige	_____
 185312 Ponto Softband, bilateral, yellow	_____
 185313 Ponto Softband, bilateral, red	_____
 185314 Ponto Softband, bilateral, pink	_____
 185315 Ponto Softband, bilateral, cerise	_____
 185316 Ponto Softband, bilateral, purple	_____
 185317 Ponto Softband, bilateral, navy	_____
 185318 Ponto Softband, bilateral, indigo	_____
 185319 Ponto Softband, bilateral, blue	_____
 185320 Ponto Softband, bilateral, turquoise	_____
 185321 Ponto Softband, bilateral, green	_____
 185322 Ponto Softband, bilateral, grey	_____
 185323 Ponto Softband, bilateral, brown	_____
 185324 Ponto Softband, bilateral, black	_____



173391
Ponto
SoundConnector

The Ponto SoundConnector allows users to use their Ponto sound processor with alternative headwear.

\$27 _____

All unilateral softbands are \$65 and all bilateral softbands are \$86.

Softband, unilateral



Softband, bilateral



		Qty:
	M50551 Soft band, pink	_____
	M50612 Soft band, blue	_____
	M50637 Soft band, black	_____
	M50638 Soft band, beige	_____
	M50639 Soft band, navy blue	_____
	M50640 Soft band, red	_____

		Qty:
	M50734 Soft band, bilateral, pink	_____
	M50735 Soft band, bilateral, blue	_____
	M50736 Soft band, bilateral, black	_____
	M50737 Soft band, bilateral, beige	_____
	M50738 Soft band, bilateral, navy blue	_____
	M50739 Soft band, bilateral, red	_____

Headband



M50308
Headband

The head band can be used by patients who are unsuitable for implantation but who can benefit from using the Ponto sound processors as a traditional bone conductor.

\$48

Because sound matters

Oticon Medical is a global company in implantable hearing solutions, dedicated to bringing the power of sound to people at every stage of life. For more than a decade, we have made bone anchored hearing systems more accessible by simplifying the treatment for physicians, audiologists, and patients alike.

We believe that patients and hearing care professionals should be able to choose the best possible solution at any time along the patient journey. We call it “Freedom of Choice” and it has always been paramount to Oticon Medical. This is the reason why our solutions are designed to be compatible whenever possible. As a result, an implant from Oticon Medical stands as a true testament to our unwavering lifelong support.

We work collaboratively to ensure that every solution we create is designed with our users’ needs in mind. We have a strong passion to provide innovative solutions and support that enhance quality of life and help people live life to the fullest – now and in the future.

Because we know how much sound matters.



 Oticon Medical LLC 580
Howard Avenue Somerset, NJ
08873 USA PHONE: +1 (888)
277-8014 FAX: +1 (732) 868-
6949
info@oticonmedicalusa.com

0825-2PS0 US